

NR3C2 Antibody



PACO10923

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

NR3C2

Uniprot

P08235

Synonyms:

nuclear receptor subfamily 3, group C, member 2; NR3C2; MCR; MGC133092; MLR; MR ;

Immunogen:

Human NR3C2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A