

NR6A1 Antibody



PACO10927

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NR6A1

Uniprot

Q15406

Synonyms:

nuclear receptor subfamily 6, group A, member 1; NR6A1; GCNF; GCNF1; NR61; RTR ;

Immunogen:

Human NR6A1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A