

NRBP2 Antibody



PACO10928

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NRBP2

Uniprot

Q9NSY0

Synonyms:

nuclear receptor binding protein 2; NRBP2; DKFZp434P086; MGC138699; TRG16; pp9320 ;

Immunogen:

Human NRBP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A