

NRD1 Antibody



PACO10930

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NRDC

Uniprot

O43847

Synonyms:

nardilysin (N-arginine dibasic convertase); NRD1; hNRD1; hNRD2 ;

Immunogen:

Human NRD1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A