

NRG1 Antibody



PACO10932

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NRG1

Uniprot

Q02297

Synonyms:

neuregulin 1; NRG1; ARIA; GGF; GGF2; HGL; HRG; HRG1; HRGA; NDF; SMDF ;

Immunogen:

Human NRG1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A