

NRG4 Antibody



PACO10933

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NRG4

Uniprot

Q8WWG1

Synonyms:

neuregulin 4; NRG4; DKFZp779N0541; DKFZp779N1944; HRG4 ;

Immunogen:

Human NRG4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A