

NRIP2 Antibody



PACO10935

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NRIP2

Uniprot

Q9BQI9

Synonyms:

nuclear receptor interacting protein 2; NRIP2; DKFZp761G1913 ;

Immunogen:

Human NRIP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A