

NSMCE2 Antibody



PACO10939

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NSMCE2

Uniprot

Q96MF7

Synonyms:

non-SMC element 2, MMS21 homolog (S. cerevisiae); NSMCE2; C8orf36; FLJ32440; MMS21; NSE2 ;

Immunogen:

Human NSMCE2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A