

PACO10941

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

NSUN5

Uniprot

Q96P11

Synonyms:

NOL1/NOP2/Sun domain family, member 5; NSUN5; FLJ10267; MGC986; NOL1; NOL1R; NSUN5A; WBSCR20; WBSCR20A; Ynl022cL; p120 ;

Immunogen:

Human NSUN5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A