

NTF4 Antibody



PACO10943

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NTF4

Uniprot

P34130

Synonyms:

neurotrophin 4; NTF4; NT-4/5; NT4; NT5; NTF5 ;

Immunogen:

Human NTF4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A