

# NUB1 Antibody



PACO10948

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NUB1

**Uniprot**

Q9Y5A7

**Synonyms:**

negative regulator of ubiquitin-like proteins 1; NUB1; BS4; NUB1L; NYREN18 ;

**Immunogen:**

Human NUB1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A