

NUBP2 Antibody



PACO10949

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NUBP2

Uniprot

Q9Y5Y2

Synonyms:

nucleotide binding protein 2 (MinD homolog, E. coli); NUBP2; CFD1; NUBP1 ;

Immunogen:

Human NUBP2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A