

NUCB1 Antibody



PACO10950

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NUCB1

Uniprot

Q02818

Synonyms:

nucleobindin 1; NUCB1; DKFZp686A15286; FLJ40471; NUC ;

Immunogen:

Human NUCB1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A