

NUF2 Antibody



PACO10962

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NUF2

Uniprot

Q9BZD4

Synonyms:

NUF2, NDC80 kinetochore complex component, homolog (S. cerevisiae); NUF2; CDCA1; NUF2R ;

Immunogen:

Human NUF2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A