

NUFIP1 Antibody



PACO10963

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NUFIP1

Uniprot

Q9UHK0

Synonyms:

nuclear fragile X mental retardation protein interacting protein 1; NUFIP1; NUFIP; bA540M5.1 ;

Immunogen:

Human NUFIP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A