

NUP50 Antibody



PACO10971

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

NUP50

Uniprot

Q9UKX7

Synonyms:

nucleoporin 50kDa; NUP50; MGC39961; NPAP60; NPAP60L ;

Immunogen:

Human NUP50.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A