

NUPL2 Antibody



PACO10977

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NUPL2

Uniprot

O15504

Synonyms:

nucleoporin like 2; NUPL2; CG1; NLP-1; NLP_1; hCG1 ;

Immunogen:

Human NUPL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A