

NUS1 Antibody



PACO10978

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NUS1

Uniprot

Q96E22

Synonyms:

nuclear undecaprenyl pyrophosphate synthase 1 homolog (S. cerevisiae); NUS1; C6orf68; MGC117249; MGC7199; MGC:7199; NgBR ;

Immunogen:

Human NUS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A