

NXF2 Antibody



PACO10983

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NXF2

Uniprot

Q9GZY0

Synonyms:

nuclear RNA export factor 2; NXF2; FLJ20416; TAPL-2 ;

Immunogen:

Human NXF2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A