

# NXF3 Antibody



PACO10984

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NXF3

**Uniprot**

Q9H4D5

**Synonyms:**

nuclear RNA export factor 3;

**Immunogen:**

Human NXF3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A