

NABP1 Antibody



PACO10987

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NABP1

Uniprot

Q96AH0

Synonyms:

oligonucleotide/oligosaccharide-binding fold containing 2A; OBFC2A; DKFZp667M1322; FLJ13624; FLJ22833; MGC111163; SSB2 ;

Immunogen:

Human OBFC2A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A