

# TENM1 Antibody



PACO10994

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TENM1

**Uniprot**

Q9UKZ4

**Synonyms:**

odz, odd Oz/ten-m homolog 1(Drosophila); ODZ1; ODZ3; TEN-M1; TNM; TNM1 ;

**Immunogen:**

Human ODZ1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A