

OGFOD2 Antibody



PACO10997

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

OGFOD2

Uniprot

Q6N063

Synonyms:

2-oxoglutarate and iron-dependent oxygenase domain containing 2; OGFOD2; DKFZp686H15154; FLJ13491; FLJ37501; MGC120434; MGC120436 ;

Immunogen:

Human OGFOD2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A