

OLFM1 Antibody



PACO11003

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

OLFM1

Uniprot

Q99784

Synonyms:

olfactomedin 1; OLFM1; AMY; NOE1; NOELIN1; OlfA ;

Immunogen:

Human OLFM1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A