

ORAI1 Antibody



PACO11011

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat, Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ORAI1

Uniprot

Q96D31

Synonyms:

ORAI calcium release-activated calcium modulator 1; ORAI1; CRACM1; FLJ14466; ORAT1; TMEM142A ;

Immunogen:

Human ORAI1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A