

PACO11014

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ORC5

Uniprot

O43913

Synonyms:

origin recognition complex, subunit 5-like (yeast); ORC5L; ORC5; ORC5P; ORC5T ;

Immunogen:

Human ORC5L

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A