

PACO11017

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat, Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

OS9

**Uniprot**

Q13438

**Synonyms:**

osteosarcoma amplified 9, endoplasmic reticulum associated protein; OS9; OS-9 ;

**Immunogen:**

Human OS9.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A