

# OSBPL11 Antibody



PACO11020

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

OSBPL11

**Uniprot**

Q9BXB4

**Synonyms:**

oxysterol binding protein-like 11; OSBPL11; FLJ13012; FLJ13164; ORP-11; ORP11; OSBP12; TCCCIA00292 ;

**Immunogen:**

Human OSBPL11.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A