

OSBPL2 Antibody



PACO11021

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

OSBPL2

Uniprot

Q9H1P3

Synonyms:

oxysterol binding protein-like 2; OSBPL2; FLJ20223; KIAA0772; MGC4307; MGC8342; ORP-2; ORP2 ;

Immunogen:

Human OSBPL2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A