

OSGEP Antibody



PACO11024

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

OSGEP

Uniprot

Q9NPF4

Synonyms:

O-sialoglycoprotein endopeptidase; OSGEP; FLJ20411; GCPL1; KAE1; OSGEP1; PRSMG1 ;

Immunogen:

Human OSGEP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A