

OSR1 Antibody



PACO11026

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

OSR1

Uniprot

Q8TAX0

Synonyms:

odd-skipped related 1 (Drosophila); OSR1; ODD ;

Immunogen:

Human OSR1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A