

OSTM1 Antibody



PACO11028

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

OSTM1

Uniprot

Q86WC4

Synonyms:

osteopetrosis associated transmembrane protein 1; OSTM1; GIPN; GL; HSPC019; OPTB5 ;

Immunogen:

Human OSTM1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A