

OTUD7B Antibody



PACO11030

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

OTUD7B

Uniprot

Q6GQQ9

Synonyms:

OTU domain containing 7B; OTUD7B; CEZANNE; ZA20D1 ;

Immunogen:

Human OTUD7B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A