

PACO11035

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

OXSM

Uniprot

Q9NWU1

Synonyms:

3-oxoacyl-ACP synthase, mitochondrial; OXSM; FASN2D; FLJ20604; KASI; KS ;

Immunogen:

Human OXSM.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A