

OXSR1 Antibody



PACO11036

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

OXSR1

Uniprot

O95747

Synonyms:

oxidative-stress responsive 1; OXSR1; KIAA1101; OSR1 ;

Immunogen:

Human OXSR1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A