

P2RX5 Antibody



PACO11038

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

P2RX5

Uniprot

Q93086

Synonyms:

purinergic receptor P2X, ligand-gated ion channel, 5; P2RX5; MGC47755; P2X5; P2X5R ;

Immunogen:

Human P2RX5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A