

# P2RY6 Antibody



PACO11041

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

P2RY6

**Uniprot**

Q15077

**Synonyms:**

pyrimidinergic receptor P2Y, G-protein coupled, 6; P2RY6; MGC15335; P2Y6 ;

**Immunogen:**

Human P2RY6.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A