

PACSIN1 Antibody



PACO11047

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PACSIN1

Uniprot

Q9BY11

Synonyms:

protein kinase C and casein kinase substrate in neurons 1; PACSIN1; KIAA1379; SDPI ;

Immunogen:

Human PACSIN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A