

PACSIN3 Antibody



PACO11049

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PACSIN3

Uniprot

Q9UKS6

Synonyms:

protein kinase C and casein kinase substrate in neurons 3; PACSIN3; SDPIII ;

Immunogen:

Human PACSIN3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A