

PAFAH1B2 Antibody



PACO11051

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

PAFAH1B2

Uniprot

P68402

Synonyms:

platelet-activating factor acetylhydrolase, isoform Ib, beta subunit 30kDa;

Immunogen:

Human PAFAH1B2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A