

PAIP1 Antibody



PACO11055

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

PAIP1

Uniprot

Q9H074

Synonyms:

poly(A) binding protein interacting protein 1; PAIP1; MGC12360 ;

Immunogen:

Human PAIP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A