

PAK1IP1 Antibody



PACO11057

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PAK1IP1

Uniprot

Q9NWT1

Synonyms:

PAK1 interacting protein 1; PAK1IP1; FLJ20624; MAK11; PIP1; RP11-421M1.5; WDR84; bA421M1.5; hPIP1 ;

Immunogen:

Human PAK1IP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A