

PANK2 Antibody



PACO11065

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PANK2

Uniprot

Q9BZ23

Synonyms:

pantothenate kinase 2; PANK2; C20orf48; FLJ17232; HARP; HSS; MGC15053; NBIA1; PKAN ;

Immunogen:

Human PANK2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A