

PANK3 Antibody



PACO11066

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PANK3

Uniprot

Q9H999

Synonyms:

pantothenate kinase 3; PANK3; FLJ12899; MGC16863 ;

Immunogen:

Human PANK3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A