

PAPOLB Antibody



PACO11068

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PAPOLB

Uniprot

Q9NRJ5

Synonyms:

poly(A) polymerase beta (testis specific); PAPOLB; PAPT; TPAP ;

Immunogen:

Human PAPOLB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A