

PAQR4 Antibody



PACO11069

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PAQR4

Uniprot

Q8N4S7

Synonyms:

progesterin and adipoQ receptor family member IV; PAQR4; FLJ30002 ;

Immunogen:

Human PAQR4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

Product Images

N/A

N/A