

PARN Antibody



PACO11072

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

PARN

Uniprot

O95453

Synonyms:

poly(A)-specific ribonuclease (deadenylation nuclease); PARN; DAN ;

Immunogen:

Human PARN.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A