

PARP15 Antibody



PACO11075

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PARP15

Uniprot

Q460N3

Synonyms:

poly (ADP-ribose) polymerase family, member 15; PARP15; BAL3; FLJ40196; FLJ40597; MGC126750; MGC126752 ;

Immunogen:

Human PARP15.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A