

PARP2 Antibody



PACO11076

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

PARP2

Uniprot

Q9UGN5

Synonyms:

poly (ADP-ribose) polymerase 2; PARP2; ADPRT2; ADPRTL2; ADPRTL3; PARP-2; pADPRT-2 ;

Immunogen:

Human PARP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A