

# PARP9 Antibody



PACO11079

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PARP9

**Uniprot**

Q8IXQ6

**Synonyms:**

poly (ADP-ribose) polymerase family, member 9; PARP9; BAL; BAL1; DKFZp666B0810; DKFZp686M15238; FLJ26637; FLJ35310; FLJ41418; FLJ43593; MGC:7868 ;

**Immunogen:**

Human PARP9.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A